



REVISIONS

REV	DR	DESCRIPTION	DATE	APPRVD
1	DAK	NEW DRAWING	9/30/04	---



NOTES:

1. DIMENSIONS ARE IN MILLIMETERS.
2. DIMENSIONS IN PARENTHESIS ARE IN INCHES.
3. REFER TO PRODUCT DESCRIPTION FOR FURTHER PRODUCT INFORMATION.
4. PROTEIN WELL VOLUME: 2uL.
5. RESERVOIR WELL VOLUME: 130uL.



CORNING
Life Sciences
45 Nagog Park
Acton, MA 01720

TOLERANCES - UNLESS SPECIFIED DIMENSIONS ARE NOMINAL FOR REFERENCE ONLY.

DRAWN DAK DATE 9/30/04
CHKD ---
MATL SEE PRODUCT DESCRIPTION

PLATE, 96 WELL, CRYSTALLOGRAPHY
SITTING DROP
3 WELLS, 1 RESERVOIR

CUSTOMER SPEC.

PROD CODE 144 SHEET 1 OF 1

DWG01859

REV 1